

SYNOLITE™ 8388-P-1

CHEMICAL/PHYSICAL NATURE

Synolite ™ 8388-P-1 is a low viscous, pre-accelerated, promoted, thixotropic, medium reactive DCPD polyester resin. An LSE system has been incorporated. Synolite™ 8388-P-1 does not contain monomeric amines.

MAJOR APPLICATIONS

Synolite™ 8388-P-1 has been especially developed for hand lay-up and spray-up applications to produce high performance constructional laminates.

PRINCIPAL PROPERTIES

Synolite™ 8388-P-1 combines good mechanical properties with optimal process ability. This resin shows very good fiber wetting and impregnation properties. Thick laminates can be made in one go due to the low exothermic heat development, good through cure and relatively low shrinkage. An optimized LSE system has been developed to get excellent Low Styrene Emission. A color change system is present for the benefit of spray-up applications. Synolite™ 8388-P-1 can be used with powder and emulsion bound mats.

APPROVALS

Synolite™ 8388-P-1 is found to comply with Det Norske Veritas' (DNV) Tentative Rules for Classification and Classification of Boats 1997, Grade 2; Det Norske Veritas' Rules for Classification of High Speed and Light Craft, Grade 2. Certificate No.: K-1720.

PRODUCT SPECIFICATIONS UPON DELIVERY


Property	Range	Unit	ТМ
Appearance	hazy	-	2265
Acid value, as such	16.5 – 23.5	mg KOH/g	2401
Viscosity, Physica, 2 s-1, 23°C	900 – 1200	mPa.s	2313
Viscosity, Physica, 20 s-1, 23°C	320 - 390	mPa.s	2313
Viscosity, Physica, 250 s-1, 23°C	200 - 230	mPa.s	2313
Solid content, IR	61 - 63	%	2033
Gel time from 25 to 35°C	22 - 26	minutes	2625
Cure time from 25°C to peak	40 - 50	minutes	2625
Peak temperature	95 - 125	°C	2625

REMARKS

Viscosity measurement: TM 2313: spindle Z2, 23°C. Curing conditions at 25°C, TM 2625: 2.0 g Butanox M-50 100 g resin




GRAPH SHOWING DYNAMIC STYRENE EMISSION



DYNAMIC STYRENE EMISSION TEST

This test is based on a preliminary test method developed to measure dynamic styrene emission. Emission is measured by an airflow of 0.4 m/s at 20°C while the surface of the resin is continuously disturbed. The description of this test method is available on request.

GRAPH SHOWING REDUCED STATIC STYRENE EMISSION




REMARKS

Measured in a stove at 23°C with certain ventilation, resins not initiated. The LSE system in Synolite™ 8388-P-1 is responsible for good low styrene emission properties.
(1) = Synolite™ 0528-P-1



GRAPH SHOWING BARCOL HARDNESS DEVELOPMENT



THROUGH CURE

Due to the special accelerator system, the final through cure of the resin in thick and thin laminates is good. In both cases the Barcol hardness after 24 hours reaches the leaflets value. Of course, the cure speed is related to the added amount of peroxide and the laminate thickness; this also defines the trimming time within broad ranges. Glass mat CSM 450 g/m^2 .

REMARKS ON CURING AGENTS

Butanox M-50 is an AKZO Nobel product, methyl ethyl ketone peroxide (MEKP).

PROCESSING

Synolite™ 8388-P-1 contains barrier forming agents to reduce emission of styrene. These agents may reduce the bonding strength of over laminates. Good strength can be obtained with over laminating of the base laminate after delayed lay-up if the surface is not too resin rich. In other cases the surface might need sanding.

GUIDELINES BEFORE USE

The resin should be conditioned at 15°C minimum before use to obtain a sufficient cure when MEKP is used as a curing system. Stir the resin mechanically before use. Do not use air bubbles to stir the resin.

STORAGE GUIDELINES

The resin should be stored indoors in the original, unopened and undamaged packaging, in a dry place at temperatures between 5°C and 30°C and the properties might change during storage. Shelf life is reduced at higher temperatures and the properties of the resin might change during storage. The shelf life of styrene containing unsaturated polyesters will be significantly reduced when exposed to light. Store in dark and in 100% light tight containers only.

MATERIAL SAFETY

A Material Safety Data Sheet of this product is available on request.

TEST METHODS

Test methods (TM) referred to in the table(s) are available on request.

The user is held to check the quality, safety and other properties of the product referred to herein. The information and recommendations in this document are to the best of our knowledge and reliable. However, no rights whatsoever can be derived from this document or the information contained therein by any party, other than those expressly accepted by a selling entity of the Aliancys group of companies ("Aliancys selling entity") in a binding sale and purchase agreement for product referred to herein. For the avoidance of doubt the Aliancys group of companies makes no warranty of any kind, express or implied, including those of merchantability and fitness for purpose. Unless explicitly agreed to otherwise in writing by the Aliancys selling entity, all offers, quotations, sales and deliveries of Aliancys of Sale of such Aliancys selling entity. Atlac®, Beyone™, Daron®, Neomould®, Neoxil®, Palatal®, Palapreg®, Synolite™, Aliancys™, the Aliancys™ logo, and the LET'S TALK/ logo are registered trademarks of Aliancys AG.